

National Institute for Public Health  
and the Environment  
*Ministry of Health, Welfare and Sport*



# Non-animal testing in the assessment of Skin sensitization

## A Sequential Testing Strategy

J.van der Veen, Emiel Rorije,  
R.Emter, A.Natsch, H.van Loveren,  
Janine Ezendam.

Dutch National Institute for Public  
Health and the Environment (RIVM)  
Bilthoven, The Netherlands

RIVM Skin sensitization testing | UK IVTS | 11-11-2014






## A Sequential Testing Strategy

- Using three “tiers” or AOP key events, combining
  - *in silico* ((Q)SAR)
  - *in chemico* (peptide reactivity)
  - *in vitro* (keratinocytes, dendritic cells)
- Stepwise approach
- Decision making involves interim decision points
- Need for additional testing is determined after each step
- Outcome:       Sensitizer / No Sensitizer
  
- Human Hazard identification





## Adverse Outcome Pathway - Allergic Contact Dermatitis





## STS – RIVM approach

- Dataset: 27 skin sensitizers; 14 non-sensitizers (incl. 2 LLNA false negatives, 4 LLNA false positives, 10 skin irritants)

- Bayesian QSAR battery (4 models)
- Direct Peptide Reactivity Assay (DPRA)

} MIE:  
peptide binding

- Keratinosens assay
- RIVM keratinocyte gene signature

} Key event  
Events in KC


- H-Clat (Dendritic cells) – published data only

→ Key event  
Events in DC

- **From simple to complex**




# HaCaT gene signature






## DNA Microarray in human keratinocytes (HaCaT)





# DNA Microarray in human keratinocytes (HaCaT)




**Predictive gene signature**

Heme oxygenase 1	Oxidative stress/Nrf2
Sulfiredoxin 1	
Stanniocalcin 2	
Adrenomedullin	
FOS	Stress response & TLR signaling
FOSL1	
Ankyrin repeat domain 37	Unknown
DNA methyltransferase 3b	DNA folding
RNA binding motif protein 5	Alternative splicing
Cyclin-dependent kinase 12	



# Predictivity gene signature

## qPCR








## TIER 1: protein binding


- Start with non-testing information
- Bayesian QSAR battery (DEREK, MultiCASE, CAESAR, OECD Toolbox)
- Equivocal results → perform DPRA





## Tier 2: Keratinocyte assays


- KeratinoSens or HaCaT gene signature
- Test **sensitizers** from tier 1 in assay with **highest PPV**
  - i.e. lowest number of false positive results
- Test **non-sensitizers** from tier 1 in assay with **highest NPV**
  - i.e. lowest number of false negative results






## Tier 3: Dendritic cell assay

- Consistent call in tier 2 → decision: sensitizer / non-sensitizer
- Equivocal call → additional information from h-CLAT required
- Decision then based on the h-CLAT outcome (majority voting)





# Results





# Predictivity

N=41	LLNA	Sequential Test Strategy
Accuracy	82.9%	100% (95.8%)
Sensitivity	92,6%	100% (96.4%)
Specificity	64,3%	100% (95.0%)
No prediction	0	2*

\* h-CLAT data was not available for all chemicals



## Conclusions and limitations

- RIVM STS provides an accurate prediction model by using the strengths of the individual tests:
  - Only using reliable QSARs data (battery) in tier 1
  - Using the most relevant *in vitro* assay in tier 2
  - Simple approach, limits (*in vitro*) testing for “known” substances
- Predictivity based on a limited number of chemicals
  - Including 4 LLNA false positives and 2 LLNA false negatives
  - Will probably decrease for a larger dataset?
- Does not provide potency information!
- Limited metabolic coverage



Contents lists available at [ScienceDirect](#)

## Regulatory Toxicology and Pharmacology

journal homepage: [www.elsevier.com/locate/yrtph](http://www.elsevier.com/locate/yrtph)



### Evaluating the performance of integrated approaches for hazard identification of skin sensitizing chemicals



Jochem W. van der Veen<sup>b,a</sup>, Emiel Rorije<sup>b</sup>, Roger Emter<sup>c</sup>, Andreas Natsch<sup>c</sup>, Henk van Loveren<sup>b,a</sup>, Janine Ezendam<sup>b,\*</sup>

<sup>a</sup> Department of Toxicogenomics, Maastricht University, PO Box 616, NL-6200 MD Maastricht, The Netherlands

<sup>b</sup> National Institute for Public Health and the Environment (RIVM), PO Box 1, NL-3720BA Bilthoven, The Netherlands

<sup>c</sup> Givaudan Schweiz AG, Ueberlandstrasse 138, CH-8600 Dübendorf, Switzerland